



Technical Data Sheet

Formaldehyde Indicator Ink for Flexography Printing Processes

REF IFVG-01

True Indicating's Formaldehyde Ink, Product Code IFVG-01, is a flexographic ink formulation used to monitor exposure to low temperature steam and formaldehyde conditions. The ink is a water based formula which contains no heavy metals or other harmful materials and has low VOC's.

Physical Properties	
Process	Low Temperature Steam Formaldehyde
Ink Carrier	Water
Packaging	HDPE or equivalent containers
pH	8.0 – 10.0
Viscosity Zahn Cup # 2	15 - 25 seconds
Ink Appearance	Thick purple-violet liquid
Printed Indicator Color	Initial Color: Violet Signal Color: Green
Anilox Roller	200 – 100 lpi
Flexographic Line Speed	≥150 ft./min.

Intended Use	
Indicator Ink is intended for manufacture of low temperature steam formaldehyde Process and Multi-Critical Variable Indicators. The Indicators are intended for use in distinguishing processed from unprocessed goods or to indicate exposure to a sterilization process at a stated value (SV). Validation of Indicators and their performance is at the discretion of the manufacturer. The Indicator Ink has been confirmed to be ISO 11140-1:2014 Type 1 compliant on Tyvek® or Film when exposed per the conditions outlined below.	
Type 1 Process Indicator	1.0 mol Formaldehyde at 70°C for 15 minutes
Printed Indicators manufactured with ink IFVG-01 can be used at temperatures between 30° and 80°C.	



True Indicating, LLC
 946 Kane Street, Suite A
 Toledo OH 43612-1372
 T: 419-476-7119

Email: info@trueindicating.com

www.TrueIndicating.com

Rev 2-02062019

® Tyvek is a registered trademark of E.I. du Pont de Nemours and Company or affiliates.

Page 1 of 2








Technical Data Sheet

Instructions for Use

IFVG-01 is intended for flexographic printing onto paper, Tyvek® or film based substrates. Formaldehyde exposure of printed ink will transition from violet to green. The shade of initial and signal color will depend on the thickness of ink printed, the substrate onto which it is printed, and the exposure conditions.

Colors shown are examples of actual printed ink using a 165 lpi Anilox roller and Tyvek®. Signal or end point color will vary due to printing conditions and exposure variables.

Initial Color	Signal Color
	

Storage and Shelf Life	
Shelf Life	7 months from date of manufacture
	15°C to 30°C (Protect from heat sources)
	Keep away from sterilant and other chemical agents. Do not reuse ink! Discard all remaining ink. Up to 5% by weight of water or ether can be added to the ink to decrease dry time.
	20% to 80% Relative Humidity

Disposal

Dispose in accordance with all federal, state, provincial, and local regulations. This material does not meet the criteria of a hazardous waste as defined under the Resource Conservation and Recovery Act (RCRA) 40 CFR 261.



True Indicating, LLC
 946 Kane Street, Suite A
 Toledo OH 43612-1372
 T: 419-476-7119

Email: info@trueindicating.com

www.TrueIndicating.com