



EMULATING INDICATOR STRIP - LABEL For Monitoring Steam Processes ISO 11140-1 TYPE 6

True Indicating Code: CSBN-16



Product Description

True Indicating Emulating Indicator Strips contain no lead or other toxic heavy metals. The Indicator Strips monitor whether sterilization conditions were met. The Indicator Strips transition from Blue to Dark with the word “NOT” being completely concealed with effective conditions. The Indicator Strips are intended for use in steam sterilizers operating at 121°C and up to 134°C. Type 6 Emulating Indicators react to all critical variables (time, temperature and moisture) and reach an end point color indicating exposure to conditions as defined by True Indicating’s stated values.

Physical Properties

Process	Steam
Indicator Strip Dimensions	19 mm x 64 mm (0.75” x 2.5”)
Packaging	250 Strips / Package
Chemical Indicator	Initial Color: Blue Signal Color: Dark The word “NOT” will be completely concealed when proper conditions have been met.

Indications for Use

Type 6 Emulating Indicator stated values:

- 121°C, 10 minutes
- 134°C, 3.5 minutes

Instructions for Use

Use a Strip in each pack, peel pouch or tray intended for steam sterilization. The Indicator strip may also be utilized as a label by peeling the liner from the strip to expose the latex free adhesive. Place strip at a location considered most challenging for steam to reach. Process packages/items as instructed in the auto-clave/sterilizer validation or manual.

Upon exposure to steam, the Indicator will transition from a shade of Blue to Dark. The transition color may vary depending on the load configuration, length and conditions of exposure. A color transition from Blue to Dark and the word “NOT” being completely concealed provides indication of effective steam sterilization conditions. If the signal color is not achieved this may suggest ideal exposure conditions were not met. Review the exposure conditions and investigate the sterilizer malfunction.




The chemical reaction which causes the color transition is a steam specific reaction and is irreversible.





Technical Data Sheet

Performance Characteristics

Result Availability	Immediately following exposure to high temperature steam processes	
Unexposed	Ineffective Exposure	Exposed to steam 134°C, 3.5 minutes
		

Colors shown are representations of unexposed, ineffective and exposed signal colors but may vary from actual use.









Color Transition is dependent upon load configuration, cycle parameters and may vary. A color transition from the Initial color to a shade of the signal color which results in the word “NOT” being completely concealed is considered an acceptable transition. If the word “NOT” remains visible, check cycle parameters and equipment function.

Compliance

ISO 11140-1:2014 Sterilization of health care products – Chemical indicators- Part 1:General requirements for Type 6 Emulating Indicators

Storage and Shelf Life

	15°C to 30°C		Keep away from sunlight
	20% to 80% relative humidity		Keep Dry
Shelf Life	3 years from the date of manufacture. The date of manufacture is based on the day the indicating ink is applied to the substrate. The remaining shelf life upon receipt will be shorter than 3 years		Do not reuse
	Keep away from sterilants. Do not use damaged indicators or indicators which have transitioned to Dark. Do not use after expiration date.		

Disposal

Discard as general waste.

