



SAFETY INFORMATION SHEET

The modified media products are not hazardous products according to the OSHA Hazard Communication Standard, 29 CFR 1910.1200. As such, a Safety Information Sheet has been established to provide general knowledge as to health and safety based upon handling, storage and use of the product.

Configuration	Product Code/REF		
Tubes	PGM-50R, PGM-100, GGM-50R, GGM-100, GGMN-100, RGM-100		
Ampules	MA-PURPLE, MA-GREEN		
Identification of Product and Company			
Product Name:	Modified Growth Medium for use with Biological Indicators and Inoculated Carriers of <i>Bacillus atrophaeus</i> , <i>Geobacillus stearothermophilus</i> , <i>Bacillus pumilus</i> or <i>Bacillus subtilis</i> for Monitoring Dry Heat, Ethylene Oxide (EO), Formaldehyde, Hydrogen Peroxide, Radiation and Steam Processes		
Product Configuration:	Tubes and Media Ampules		
Intended Use:	Accessory to Indicator		
NFPA 704 Hazard Rating	Health: 0	Fire: 0	Reactivity: 0
Prepared By:	J. Wheeler	Email:	info@trueindicating.com
Company:	True Indicating, LLC 946 Kane St, Suite A, Toledo Ohio 43612		
	Emergency Telephone No.: 1-419-476-7119 (True Indicating); 1-800-424-9300 (CHEMTREC)		
Product Composition			
Product consists of Soybean Casein Digest Broth (SCDB), Deionized (DI) Water and a pH Indicator (less than 0.5% by volume) which transitions from the initial color to yellow when bacterial growth occurs.			
Product is free of latex and natural dry rubber products.			
First Aid Measures			
Eyes:	If in eyes, remove contact lenses where applicable and rinse with water for 15 minutes.		
Skin:	Wash with soap and water as normal hygiene dictates.		
Ingestion:	Not a likely route of exposure, however, if ingestion occurs, give water to drink and seek medical advice.		
Inhalation:	No safety concerns associated with inhalation.		
Fire Fighting Measures			
Suitable Extinguishing Media:	Use water spray, alcohol resistant foam, dry chemical or carbon dioxide.		
Fire Hazard:	Not flammable	Explosion Hazard:	Product is not explosive
Handling			
Precautions for Safe Handling:			
Handle in accordance with good industrial hygiene and safety procedures. Wash hands and other exposed areas with soap and water before eating, drinking or smoking and again when leaving work.			



SAFETY INFORMATION SHEET

Shipping Information		
Hazard Category:	Non-hazardous	
Conditions/Restrictions:	None	
Harmonized Tariff Code:		
Product REF	Code	Description
All Listed	3821.00.0000	Prepared culture media for development or maintenance of micro-organisms or plant, human or animal cells