



This lot of product meets the accepted performance criteria recommended in ISO 11140-1: 2014 as a Type 1 Process Indicator.


Product Name: Indicator for Monitoring Ethylene Oxide (EO) Processes


For use in Monitoring: Ethylene Oxide (EO)

Appearance: EO Indicator Tape, 19 mm x 50m





LOT BC210000

REF EO-TAPE

 2021-01-01

 2023-01-01

Characteristic	Specification	Result
Ethylene Oxide (100%) 54°C, 20 minutes, 600 mg/L, 60% RH	Red/Orange to Green Color Transition	Pass
Date of Product Performance Verification:	January 2021	

Storage and Shelf Life			
	Room Temperature (15° – 30°C)		Keep away from sunlight
	20% to 80% Relative Humidity		Protect from heat, radioactive sources, and sterilizing agents
Shelf Life	24 Months from the date of manufacture		

Disposal: Discard as general waste.

May 20, 2021

Quality Approval

Date