



This lot of product meets the accepted performance criteria recommended in ISO 11140-1: 2014 as a Type 1 Process Indicator.

Product Name: Indicator for Monitoring Steam Processes



B182727



CSWN-1T



2020-10-17



2018-10-17

For use in Monitoring: Steam

Appearance: Steam Indicator Tape, 18 mm x 50 m

Performance Characteristics

Characteristic	Specification	Result
Steam 121°C, 10 minutes	Light to Dark Color Transition	Pass
Date of Product Performance Verification:	October 2018	

Storage and Shelf Life

	Room Temperature (15° – 30°C)		Keep away from sunlight
	20% to 80% Relative Humidity		Protect from heat, radioactive sources, and sterilizing agents
Shelf Life	24 Months from the date of manufacture		

Disposal: Discard as general waste.

\_\_\_\_\_  
Quality Approval

\_\_\_\_\_  
Date