



This lot of product meets the accepted performance criteria recommended in ISO 11140-1: 2014 as a Type 5 Integrating Indicator.

Product Name: Indicator for Monitoring Steam Processes

LOT S180531 REF CSPN-15 31-MAY-2021

31-MAY-2018 For use in Monitoring: Steam

Appearance: Indicator Strip-Label, 19 mm x 64 mm

Performance Characteristics

Table with 3 columns: Characteristic, Specification, Result. Row 1: Steam 121°C, 16.5 minutes or 135°C, ≥1.2 minutes; Pink to Dark Color Transition Completely Covering the Word "Not"; Pass. Row 2: Date of Product Performance Verification: May 2018

Storage and Shelf Life

Table with 4 columns: Storage condition icon, Temperature/Relative Humidity, Hazard icon, Instruction. Row 1: Room Temperature (15° – 30°C), Keep away from sunlight. Row 2: 20% to 80% Relative Humidity, Protect from heat, radioactive sources, and sterilizing agents. Row 3: Shelf Life: 36 Months from the date of manufacture

Disposal: Discard as general waste.

Quality Approval

Date