




This lot of product meets the accepted performance criteria recommended in ISO 11140-1: 2014 as a Type 4 Multi-critical Process Variable Indicator.

Product Name: Indicator for Monitoring Steam Processes

**LOT** S180531                                      **REF** CSYN-14                                       31-MAY-2021





 31-MAY-2018                                      For use in Monitoring:      Steam

Appearance:      Indicator Strip-Label, 19 mm x 64 mm

PERFORMANCE CHARACTERISTICS

Characteristic	Specification	Result
Steam 121°C, 10 minutes or 134°C for 3.5 minutes	Yellow to Dark Color Transition Completely Covered the Word "Not"	Pass
Date of Product Performance Verification:	MAY 2018	

Storage and Shelf Life

	Room Temperature (15° – 30°C)		Keep away from sunlight.
	20% to 80% Relative Humidity		Protect from heat, radioactive sources, and sterilizing agents.
Shelf Life	36 Months from the date of manufacture		

Disposal: Discard as standard waste.

\_\_\_\_\_  
Quality Approval

\_\_\_\_\_  
Date