



This lot of product meets the accepted performance criteria recommended in ISO 11140-1: 2014 as a Type 1 Process Indicator.

Product Name: Chemical Process Indicator for Monitoring Dry Heat Processes


For use in Monitoring: Dry Heat

Appearance: Blue to Dark, 19 mm x 32 mm Label





LOT D210101

REF CDBN-2

 2021-01-01

 2025-01-01

Characteristic	Specification	Result
Dry Heat 160°C, 40 minutes	Blue to Dark Color Transition	Pass
Date of Product Performance Verification	January 2021	

Storage and Shelf Life			
	Room Temperature (15° – 30°C)		Keep away from sunlight
	20% to 80% Relative Humidity		Protect from heat, radioactive sources, and sterilizing agents.
Shelf Life	4 years from the date of manufacture		

Disposal: Discard as general waste.

July 20, 2021

Quality Approval

Date