



This lot of product meets the accepted performance criteria recommended in ISO 11140-1: 2014 as a Type 4 Multi-Critical Process Variable Indicator.

Product Name: Type 4 Multi-Critical Process Variable Indicator

For use in Monitoring: Vaporized Hydrogen Peroxide (VH<sub>2</sub>O<sub>2</sub>) Sterilization Processes

Appearance: Indicator Strip Label, 19 mm x 64 mm

**LOT** P200624

**REF** CPHY-14

24-JUN-2020

24-JUN-2023

Characteristic	Specification	Result
Vaporized Hydrogen Peroxide (VH <sub>2</sub> O <sub>2</sub> ) 50°C, 6 minutes, 2.3 mg/L	Red to Yellow Color Transition	Pass
Date of Product Performance Verification:	June 2020	

Storage and Shelf Life			
	Room Temperature (15° – 30°C)		Keep away from sunlight
	20% to 80% Relative Humidity		Protect from heat, radioactive sources, and sterilizing agents
Shelf Life	36 Months from the date of manufacture		

Disposal: Discard as general waste.

\_\_\_\_\_  
Quality Approval

\_\_\_\_\_  
Date