



This lot of product meets the accepted performance criteria recommended in ISO 11140-1: 2014 as a Type 1 Process Indicator.

Product Name: Chemical Process Indicator for Monitoring Vaporized Hydrogen Peroxide (VH<sub>2</sub>O<sub>2</sub>) Processes

For use in Monitoring: Vaporized Hydrogen Peroxide (VH<sub>2</sub>O<sub>2</sub>)

Appearance: Yellow to Blue, 12 mm Imprinted Dot

**LOT** P220222

**REF** CPYB-1

2022-02-22

2026-02-22

Characteristic	Specification	Result
VH <sub>2</sub> O <sub>2</sub> 50°C, 2.3 mg/L VH <sub>2</sub> O <sub>2</sub> , 6 minutes	Yellow to Blue Color Transition	Pass
Date of Product Performance Verification:	February 2022	

Storage and Shelf Life			
	Room Temperature (15° – 30°C)		Keep away from sunlight
	20% to 80% Relative Humidity		Protect from heat, radioactive sources, and sterilizing agents
Shelf Life	4 years from the date of manufacture.		

Disposal: Discard as general waste.

August 2, 2022

\_\_\_\_\_  
Quality Approval

\_\_\_\_\_  
Date