



SAFETY INFORMATION SHEET

The Biological Indicator and Inoculated Carrier products are not hazardous products according to the OSHA Hazard Communication Standard, 29 CFR 1910.1200. As such, a Safety Information Sheet has been established to provide general knowledge as to health and safety based upon handling, storage and use of the product.

Configuration		Product Code/REF	
Spore Discs	DA-06, DA3-06, DT-05, DT-06, DT3-06, SDAM-06, SDAT-06, SDGT-0x, SDTM-06, SDTT-0x		
Mini-Spore Strips	MAB-06, MA-06, MTB-06, MT-06		
Spore Strips	SA-0x, SB5-0x, SBP-0x, ST-0x		
Spore Threads	TAB-06, TA-06, TTB-06, TT-06		
Spore Wires	WAB-06, WA-06, WTB-06, WT-06, WTH-06		
Identification of Product and Company			
Product Name:	Biological Indicators and Inoculated Carriers of <i>Bacillus atrophaeus</i> , <i>Geobacillus stearothermophilus</i> , <i>Bacillus pumilus</i> or <i>Bacillus subtilis</i> for Monitoring Dry Heat, Ethylene Oxide (EO), Formaldehyde, Hydrogen Peroxide, Radiation and Steam Processes		
Product Configuration:	Spore Discs, Spore Strips, Spore Threads and Spore Wires		
Intended Use:	Indicator		
NFPA 704 Hazard Rating	Health: 0	Fire: 0	Reactivity: 0
Prepared By:	J. Wheeler	Email:	info@trueindicating.com
Company:	True Indicating, LLC 946 Kane St, Suite A, Toledo Ohio 43612		
	Emergency Telephone No.: 1-419-476-7119 (True Indicating) 1-800-424-9300 (CHEMTREC)		
Product Composition			
Product consists of filter paper or a steel carrier impregnated with <i>Bacillus atrophaeus</i> , <i>Geobacillus stearothermophilus</i> , <i>Bacillus pumilus</i> or <i>Bacillus subtilis</i> spores in known concentrations.			
Product is free of latex and natural dry rubber products.			
First Aid Measures			
Eyes:	Not applicable		
Skin:	No safety concerns associated with skin contact; wash with soap and water as normal hygiene dictates.		
Ingestion:	Not a likely route of exposure, however, if ingestion occurs, give water to drink and seek medical advice.		
Inhalation:	No safety concerns associated with inhalation.		
Fire Fighting Measures			
Suitable Extinguishing Media:	Use water spray, alcohol resistant foam, dry chemical or carbon dioxide		
Fire Hazard:	Not flammable	Explosion Hazard:	Not explosive



SAFETY INFORMATION SHEET

Handling	
Precautions for Safe Handling: Handle in accordance with good industrial hygiene and safety procedures. Wash hands and other exposed areas with soap and water before eating, drinking or smoking and again when leaving work.	
Methods for Containment and Cleaning Up	
For Containment:	Not Applicable
Methods for Cleaning Up:	Clean up spills immediately and dispose of waste safely.
Shipping Information	
Hazard Category:	Non-hazardous
Conditions/Restrictions:	None
Harmonized Tariff Code:	3002.90.5150