



SAFETY INFORMATION SHEET

The Instant 20s products are not hazardous products according to the OSHA Hazard Communication Standard, 29 CFR 1910.1200. As such, a Safety Information Sheet has been established to provide general knowledge as to health and safety based upon handling, storage and use of the product.

	Product Code/REF		
Instant 20s	ISCS-05		
Identification of Product and Company			
Product Name:	Instant 20s Indicators		
Product Configuration:	Enzymatic Indicators (EI)		
Intended Use:	Instant 20s Indicator, a multiple, interactive, bacterial enzyme indicator, is used for monitoring saturated steam sterilization processes operating at: 121°C, 30 minutes (Gravity), 132°C, 4 minutes (Pre-Vacuum), 135°C, 3 minutes (Pre-Vacuum).		
NFPA 704 Hazard Rating	Health: 0	Fire: 0	Reactivity: 0
Prepared By:	J. Wheeler	Email:	info@trueindicating.com
Company:	True Indicating, LLC 946 Kane St, Suite A, Toledo Ohio 43612		
	Emergency Telephone No.: 1-419-476-7119 (True Indicating); 1-800-424-9300 (CHEMTREC)		
Product Composition			
Product is free of latex and natural dry rubber products.			
The Instant 20s Indicator consists of a polypropylene vial and cap, a silicone vented base with foam plug insert and a biological tablet produced using cellulose, and enzymes from <i>Geobacillus stearothermophilus</i> ATCC® 7953.			
First Aid Measures			
Eyes:	Not applicable		
Skin:	No safety concerns associated with skin contact; wash with soap and water as normal hygiene dictates.		
Ingestion:	Not a likely route of exposure, however, if ingestion occurs, give water to drink and seek medical advice.		
Inhalation:	No safety concerns associated with inhalation.		
Fire Fighting Measures			
Suitable Extinguishing Media:	Use water spray, alcohol resistant foam, dry chemical or carbon dioxide.		
Fire Hazard:	Not flammable	Explosion Hazard:	Product is not explosive
Handling			
Precautions for Safe Handling:			
Handle in accordance with good industrial hygiene and safety procedures. Wash hands and other exposed areas with soap and water before eating, drinking or smoking and again when leaving work.			
Shipping Information			
Hazard Category:	Non-hazardous		
Conditions/Restrictions:	None		
Harmonized Tariff Code:			
Product REF	Code	Description	
ISCS-05	3507.90.7000	Enzymes, prepared enzymes not elsewhere specified or included, Other, Other	