



This lot of product meets the accepted performance criteria recommended in ISO 11140-1: 2014 as a Type 1 Process Indicator.


Product Name: Chemical Process Indicator for Monitoring Ethylene Oxide (EO) Processes


For use in Monitoring: Ethylene Oxide (EO)

Appearance: Violet to Green, 12 mm Imprinted Dot





LOT E210101

REF CEVG-3

 01-JAN-2021

 01-JAN-2026

Characteristic	Specification	Result
Ethylene Oxide (100%) 600 mg/L, 54°C, 60% RH, 20 minutes	Violet to Green Color Transition	Pass
Date of Product Performance Verification	January 2021	

Storage and Shelf Life			
	Room Temperature (15° – 30°C)		Keep away from sunlight
	20% to 80% Relative Humidity		Protect from heat, radioactive sources, and sterilizing agents
Shelf Life	5 Years from the date of manufacture		

Disposal: Discard as general waste.

October 6, 2021

Quality Approval

Date